

VARIABLE ANGLE REFLECTION ACCESSORY WITH VW ABSOLUTE REFLECTANCE STAGE NO. 21149

Comparison of Absolute Reflectance Methods: VW vs Reflectance Relative to a Standard



Figure 1. The Variable Angle Reflection Accessory with its Variable Angle Stage.




Figure 2. The 12° VW Absolute Reflection Stage, in its V (upper) and W (lower) configurations.

INTRODUCTION

Absolute reflectance is valuable in the optics and semiconductor industries. It can be used to characterize mirrors and other optical components in addition to coatings thereon. Absolute reflection measurements are necessary to determine the refractive indices and optical constants of materials.

Several different methods can be used to measure absolute reflectance. The two most common techniques are the VW method and comparison to a standard.

For the VW measurements, the background is collected with the stage in V mode and the sample collected in W mode and then the reflectance is calculated from the measurement. The VW offers the advantage of introducing no additional optics between the sample and reference measurement, thereby minimizing errors. For the measurement relative to a standard, an additional optic – the standard – is introduced. So the measurement depends on the refractive index of the reference.

The two methods are compared here to see if they give equivalent results.

EXPERIMENTAL

All spectra were collected on a commercial FT-IR spectrometer with the Harrick [Variable Angle Reflection Accessory](#) (see Figure 1) installed, with its VW absolute reflectance stage with a 12° incident angle (see Figure 2). Note that the reflectance depends on the incident angle and polarization, as well as the properties of the sample. A 12° incident angle is considered sufficiently low to avoid a significant contribution from the incident polarization so the measurements herein were done without a polarizer.




Figure 3. Absolute Reflectance Spectrum of SiO₂ measured using the VW method (red) and relative to a reference (blue).

Two sets of measurements were recorded. For the VW measurements, the background was collected with the stage in V mode and the sample collected in W mode. The square root of the resulting spectrum was taken to extract the reflectance from the collected data. For the relative reflectance measurements, both the background and the sample were recorded with the stage in the V mode. The background was collected using a Ge reflectance reference (Harrick Scientific, P/N RRF-00J).

All spectra were signal averaged over 64 scans over the 4000 cm⁻¹ to 400 cm⁻¹ region. A background was collected before every sample spectrum. The gain and optical velocity were adjusted to maximize the energy on the TGS detector while avoiding detector saturation.

Three materials were investigated: UV-grade fused silica (SiO₂) and calcium fluoride (CaF₂).

RESULTS AND DISCUSSION

The absolute reflectance of SiO₂ and CaF₂ are shown in Figures 3 and 4 respectively. The data in red was collected using the VW method while that in blue was collected using the reflectance reference. In all cases, the two methods give virtually identical results.




Figure 4. Absolute reflectance of CaF₂ measured using the VW method (red) and relative to a reference (blue).

The spectra for both materials are very similar. This is not surprising because both materials have the same refractive index (SiO₂ and CaF₂, n=1.4). Note the slightly higher than typical noise level in the relative reflectance measurements for both materials, which may be attributed to the reference used. Using a reference with a refractive index more closely matched to that of the sample might reduce the noise.

CONCLUSION

Both the VW and the relative reflectance absolute reflectance methods give comparable results in terms of the measured reflectance.