



DIAPHRAGM OPTICAL ATTENUATOR

APPLICATIONS

- ▶ Handy for aligning FTIR spectrometers.
- ▶ Useful for precisely controlling the beam size.
- ▶ Eliminates stray radiation.

FEATURES

- ▶ Continuously variable aperture.
- ▶ Diameter of the iris can be changed from 25mm to 0.5mm.
- ▶ Scale for repeatable setting.

INCLUDES

- ▶ Diaphragm Optical Attenuator.
- ▶ Slide plate mounted.



ORDERING INFORMATION

Diaphragm Optical Attenuator.....

CATALOG No.
OAD-S1G