



UNIVERSAL TRANSMISSION HOLDER

Our Universal Transmission Holder is ideal for securing a variety of samples for analysis by transmission spectroscopy. Its sliding MagnaGrip™ plate firmly holds both thin films and samples up to 38mm thick in the beam of the FT-IR spectrometer.

APPLICATIONS

- ▶ Convenient for holding thin polymer films, windows, pellets and KBr pellet holders firmly in place for spectroscopic analysis.

FEATURES

- ▶ Highly versatile transmission sample holder.
- ▶ MagnaGrip™ plate secures thin and flexible samples against a sample slide plate.
- ▶ MagnaGrip™ plate slides on rods to accommodate samples up to 38mm long.
- ▶ 20mm round aperture.
- ▶ Compatible with standard FT-IR and UV-Vis sample slide holders.

INCLUDES

- ▶ Universal Transmission Holder.
- ▶ Inserts with 13mm, 6mm, and 3mm apertures for securing small samples.
- ▶ Sample slide plate mounting.



ORDERING INFORMATION

Universal Transmission Holder

CATALOG NO.
HUT-TS1

Harrick Scientific Products, Inc.
141 Tompkins Ave, 2nd floor, Pleasantville NY 10570
Ph: 800-248-3847, FAX: 914-727-7209, web site: www.harricksci.com, e-mail: info@harricksci.com