



## SPECULAR REFLECTANCE ACCESSORIES (SLIDE PLATE MOUNTED)

### APPLICATIONS

- ▶ Reflectance measurements of films on metallic substrates, epitaxial films, anti-reflection and enhanced reflection coatings, etc.
- ▶ Reflectance measurements of opaque materials.

### FEATURES

- ▶ Convenient horizontal sampling surface.
- ▶ High throughput.
- ▶ Ease of use.
- ▶ Interchangeable masks for examining small samples or small areas on large samples.
- ▶ Convenient slide plate mounting for multi-use spectrometers.
- ▶ Fixed 12.5°, 22.5°, 30°, 45° or 60° angles of incidence.
- ▶ Near-normal (8°) floor mounted model also available.



### INCLUDES

- ▶ Specular reflectance accessory.
- ▶ Mask set with masks of three different sizes: 3, 6 and 13 mm in diameter.
- ▶ Alignment mirror.

### ORDERING INFORMATION

12.5° Specular Reflectance Accessory.....	ERA-12G
22.5° Specular Reflectance Accessory.....	ERA-22G
30° Specular Reflectance Accessory.....	ERA-30G
45° Specular Reflectance Accessory.....	ERA-45G
60° Specular Reflectance Accessory.....	ERA-60G
Near-Normal (8°) Specular Reflectance Accessory (floor mounted).....	ERA-8-XXX

*XXX denotes the spectrometer code*

### OPTIONAL ACCESSORIES FOR ALL MODELS EXCEPT NEAR-NORMAL

Glan-Taylor Polarizer Mount.....	ERA-PGT-HDR
Glan-Taylor Polarizer.....	PGT-S1V
Depolarizer and Wire Grid Polarizer (PWG-U1R) Mount.....	ERA-PWG-HDR
Depolarizer.....	DPS-R4V
Wire Grid Polarizer.....	PWG-U1R