



OPTICAL POLISHING KIT

Useful for grinding, buffing and polishing many common optical materials that do not require pitch polishing. Optical polishing is an art; it takes expertise to obtain high quality polish.

APPLICATIONS

- ▶ A convenient method for grinding, polishing, and buffing NaCl, KBr, ZnSe, and KRS-5 windows and ATR plates.

FEATURES

- ▶ Contains all the items needed to improve the surface finish of optical windows.

INCLUDES

- ▶ Two ground glass plates.
- ▶ Two self-adhering polishing pads.
- ▶ #400 coarse grinding compound.
- ▶ #600 fine grinding compound.
- ▶ Cerium oxide polishing compound.
- ▶ KRS-5 polishing compound.
- ▶ Finger cots.
- ▶ One camel hair brush.
- ▶ 2 oz. empty bottle.
- ▶ Five pouring spouts.



ORDERING INFORMATION

Optical Polishing Kit

CATALOG NO.
OPK-1XX

OPTIONS AND REPLACEMENT PARTS

Ground Glass Plate	OPK-GGP
Polishing Pad (self-adhering)	OPK-MPP
#400 Coarse Grinding Compound	OPK-CGC
#600 Fine Grinding Compound	OPK-FGC
Cerium Oxide Polishing Compound	OPK-REP
Brush	OPK-CHB
Polishing Compound for KRS-5	OPK-KPC
Aluminum Oxide Polishing Compound	OPK-APC
Finger Cots (100 pcs)	KIT-COT