



OPTICAL TABLE SLIDE PLATE HOLDER

APPLICATIONS

- ▶ Convenient for holding sample slide plate mounted accessories.
- ▶ Firmly holds Demountable Liquid Cells, Demountable Gas Cells, sample holders, attenuators, polarizers and depolarizers onto optical tables, benches and breadboards.

FEATURES

- ▶ Mounts standard 2" x 3" x 0.125" slide plate mounted sampling accessories.
- ▶ Compatible with the standard posts used with optical tables.

INCLUDES

- ▶ Sample Slide Holder with three Metric M6 or English 1/4-20 mounting holes.
- ▶ Screw for securing the slide plate.



ORDERING INFORMATION

Sample Slide Holder, English
Sample Slide Holder, Metric

CATALOG NO.
HSS-OTE
HSS-OTM