



RAIL MOUNT

APPLICATIONS


- ▶ System for mounting accessories in sample compartments securely and precisely.
- ▶ Simplifies operations where other items, such as polarizers, must be placed co-linear with the sample beam.

FEATURES

- ▶ Enables up to three slide plate holders to be installed in the beam path at one time.
- ▶ Permits translation of accessories along the beam path to optimize performance.
- ▶ Compatible and included with some Harrick Scientific accessories.
- ▶ Available for most FTIR and UV-Vis spectrometers.

INCLUDES

- ▶ Rail Mount.
- ▶ Mating hardware for the specified spectrometer.



ORDERING INFORMATION

Rail Mount.....

CATALOG NO.
HRM-XXX

OPTIONS AND REPLACEMENT PARTS

Sample Slide Holder.....

HSS-XXX