



SAMPLE SLIDE HOLDER

APPLICATIONS

- ▶ Convenient for holding sample slide plate mounted accessories.
- ▶ Firmly holds Demountable Liquid Cells, Demountable Gas Cells, sample holders, attenuators, polarizers and depolarizers into FTIR and UV-Vis spectrometers.

FEATURES

- ▶ Mounts standard 2" x 3" x 0.125" slide plate mounted sample holders.
- ▶ Compatible with Harrick Scientific's Rail Mount.
- ▶ Available for FTIR and UV-Vis spectrometers.

INCLUDES

- ▶ Sample Slide Holder.
- ▶ Screws for securing the slide plate.



ORDERING INFORMATION

Sample Slide Holder.....

CATALOG NO.
HSS-XXX

OPTIONS AND REPLACEMENT PARTS

Rail Mount.....

HRM-XXX