




## NEAR-NORMAL INCIDENCE REFLECTANCE ACCESSORY

### APPLICATIONS

- ▶ Reflectance measurements of films on metallic substrates, epitaxial films, anti-reflection and enhanced reflection coatings.
- ▶ Suitable for use with opaque materials.

### FEATURES

- ▶ Fixed 8° angle of incidence.
- ▶ Horizontal sampling surface.
- ▶ Repeatable 3-point sample support.
- ▶ High throughput for optimum performance.
- ▶ Ease of use.
- ▶ Interchangeable masks for examining small samples or small areas on large samples.
- ▶ Slide plate mounted models available for 12.5°, 22.5°, 30°, 45° and 60° angles of incidence.



### INCLUDES

- ▶ Specular reflection accessory.
- ▶ Mask set with masks of three different sizes: 3, 6 and 13 mm in diameter.
- ▶ Alignment mirror.
- ▶ Hardware to install into the spectrometer.

### ORDERING INFORMATION

	CATALOG NO.
Near-Normal Incidence Reflectance Accessory.....	ERA-8-XXX
12.5° Specular Reflectance Accessory (slide-plate mounted).....	ERA-12G
22.5° Specular Reflectance Accessory (slide-plate mounted).....	ERA-22G
30° Specular Reflectance Accessory (slide-plate mounted).....	ERA-30G
45° Specular Reflectance Accessory (slide-plate mounted).....	ERA-45G
60° Specular Reflectance Accessory (slide-plate mounted).....	ERA-60G

XXX denotes the spectrometer code